

Safekeeper^{COLD}



Safekeeper COLD Evidence Refrigerator

- For Temperature Controlled Storage with Multiple Layers of Security for Pharmaceutical and Forensic Applications



Safekeeper[®] COLD



"The World's Most Extensive Selection of Containment Solutions."



CONTENTS:

- Product Overview (p.2)
- Specifications (p.3)

INTRODUCTION

When evidence must be stored below room temperature to maintain integrity, use of a household or shared refrigerator may not be satisfactory to provide uniform storage conditions or security in the chain of custody.

The Safekeeper COLD is a unique refrigerator designed to provide temperature controlled storage within a single cabinet with multiple layers of security.

APPLICATIONS

Pharmaceutical Storage \ Forensic Investigation and Processing



KEY FEATURES

Compact Footprint. The Safekeeper COLD is less than 24" wide, which allows placement in nearly any location.

Customizable Door. Reversible door with option to customize with a 304-stainless steel wrap.

100% CFC Free. Environmentally friendly design with no ozone-damaging chemicals.

PERFORMANCE ADVANTAGE

The 10.1 cu. ft. Safekeeper COLD refrigerator features nine locking interior compartments to keep evidence separated and secure. Each compartment is completely enclosed to keep contents from commingling and is protected by a combination lock with user-defined codes. Two master keys are provided in case of lost combinations or administrative overrides.

An internal blower system circulates conditioned air throughout the cabinet. Open wire shelving enhances uniformity at all storage levels. A fully automatic defrost system activates on demand to optimize performance of the refrigeration system.



CONTENTS:

Product Overview (p.2)

Specifications (p.3)


Safekeeper^{COLD}

ALL Refrigerator

SPECIFICATIONS

3

FRONT VIEW



TOP VIEW



MODEL	VOLUME	DIMENSIONS		WEIGHT (LBS/KG)	
		Internal (W x D x H)	External (W x D x H)	Net	Ship
Safekeeper COLD	10.1 / 286	20.5" x 17.75" x 48.25" 520 x 450 x 1230 mm	23.63" x 23" x 55.5" 600 x 584 x 1410 mm	140 / 64	145 / 66

Options and Accessories

- Optional 304 stainless steel door wrap.

PRODUCT SPECIFICATIONS

COOLING SYSTEM

Safekeeper COLD

Airflow <... Chamber air circulated through the wire shelves and along the door surface by an internal blower with positive airflow ...>

CONSTRUCTION

Safekeeper COLD

Finish <... Coated metal interior and exterior ...>

Door <... White metal door, right hand swing; field reversible. Optional 304 stainless steel wrap is available. ...>

Shelves <... 4 metal wire shelves, 9 lockable compartments ...>

Electrical <... 115 V AC/60 Hz ...>

Temperature Indicator <... Dial thermometer ...>

Compressor <... 115 V AC/60 Hz; 3.0 oz. capacity (R134a refrigerant required) ...>

Leveling Feet <... 2 ...>

CONTENTS:

Product Overview (p.2)

● Specifications (p.3)

Safekeeper^{COLD}

ALL Refrigerator

SPECIFICATIONS

4

WARRANTY

This product is protected by the Air Science Legacy Limited Lifetime Warranty™.



For details visit the [Warranty section](#) of our website.



120 6th Street \ Fort Myers, FL 33907
T. 239-489-0024 \ **Toll Free.** 800-306-0656 \ **F.** 800-306-0677
www.airscience.com

The information contained in this manual and the accompanying product are copyrighted and all rights are reserved by Air Science. Air Science reserves the right to make periodic minor design changes without obligation to notify any person or entity of such change.

